

DKK2-dickkopf homolog 2 (*Xenopus laevis*) Gene

Reference: PTXBC126328

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DKK-2; dickkopf-related protein 2; dickkopf 2 homolog; dickkopf homolog 2; dickkopf-2; hDkk-2; dickkopf WNT signaling pathway inhibitor 2		

For research use only.

GENE DESCRIPTION	dickkopf homolog 2 (<i>Xenopus laevis</i>)
GENE ID	27123
ORIGIN SPECIES	Human
ACCESSIONS	BC126328
VECTOR	pDONR223
NCBI SYMBOL	DKK2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccggttgatcggagcaaggattcgtctgctgctcctactggccggtgctgatgggagagctcacagatcggcagttcggggccaaactcaatccatcaagtcctctctgggctggggagacgcctggtagccccaatgatctgcggcatgtaccaagg
actggcattcggcggcagtaagaaggcgaacacctggggcaggcctacccctgtagcagtgataaggagtggaagtgggaggtattgccacagtcgccacccaaggatcatcgccctgcatgggtgtctggagaaaaagaagcgtgccaccgagatggcatgtgctgccccagt
accgctgcaataatggcatctgtatccagttactgaaagcatcttaacccctcacatccggctctggtgactcagcacagagatcgaaccacgggtcattactcaaacatgacttggatggcagaatctaggagaccacacactaagatgtcacatataaaggcatgaagg
agacccctgcctacgatcatcagactgattgaagggttttctgtgctcgtcattttggacaaaatctgcaaacagtgctccatcagggggaagtctgtaccaaaacgcaagaagggttctcatgggtggaaattttccagcgttcgactgtgcgaaggcctgtctgcaaaagta
tggaagatgccactactcctcaaaaggcagactcatgtgtgcagaaatttga