

## C7orf38-chromosome 7 open reading frame 38 Gene

Reference: PTXBC126337

|                               |                                                                        |                     |                |
|-------------------------------|------------------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | C7orf38; protein FAM200A; family with sequence similarity 200 member A |                     |                |
| <b>For research use only.</b> |                                                                        |                     |                |

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | chromosome 7 open reading frame 38                                                               |
| <b>GENE ID</b>          | 221786                                                                                           |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC126337                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | C7orf38                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgactcctgaatcaagggaactacagatttctccagggtaccaggagatggaaggcatcgtgatagtgaagggtgaggaggagatgaagaagaccattttcaaaaggaaagaaacaaagtagagtcacgccacaagttctcagtcgctctacaactatgaatgagagagcc  
ttattatcatcgtattttagtgcataatagatggcctcaaaagagaaaatggctcacacagcggctgaaaaaattatctccagcatgtatggacatggtagcacaattttgatgacaaatcagctgataaactaagaactatactcttagtgataatacaatctcgtgaatctgtacgattgca  
aaacatttggaaagcaatgcttattacacggctcagtcggtatagactttgcaatcacaactcgtatgagagcactgatatgcaagttgcccacactcttggtttatgtcagatattgtggcaagatgattttagaggatctctatgttttaaaattaaattcacatataactggattgatttatt  
tactgaattagaaaactgcctcttggtcagtatataaattaaactggaacacttgaaggaatttcaagtgatggaacgcaaatatgaccggaacacacagcagacttactgaaaaattgttagaagcaaccacacaatgctgttggaaatcactgttttattcatcagaagccttggatcca  
aagaaattccacaaagctgatggatgattgaaaaatgcagtgaacactgttaattttataaaggaagtcactgaatagccgacttctcgaatattttgtcagagattggagtgaaccacaccacttattgtttcatacagaagtcgttgccttctcaaggaaaaatgattgacagagatat  
gaactcagggaacgagatttatttctcgttgaagaaatcctcatttggcaaatattttgaagacgacatttggtaacaaaattggcatatttaagtgtatattttggcattcttaataaagcctgaaatgcaggggaaaaacaaatgatatttctgattctgaacatattcttaggattcc  
aaaagacgttatttatttgcgaagcaagacttaaaagtaaccgcccctagctactatattttccaacattattgcaacacatcgaagagaacatttaataagagactgcttaaaagaaataaattagagatattgttgcattcactcttcttctcaaaacttttaattattcttccggaagagaatt  
tgaatcattaaagaaaatatttggatgaagatccatttgcctttcaaaaccagaatcaataattgagttaaacttgagcctgaagagaagatgaattattgcagctcagttcatcattacactaaagaattattataagatttaagttatcagcatttggattaaagattaaagatgactttcca  
ctgctaagtaggaagagatatttctgttactacattcacaactacatatttctgtgaactagatttcaatcttgacacgggttaaaacaaagaagagaataagctcaatattgcaccagatatgcgggtagcattatcttcatgttctcactggaaggaaactatgaacagacaagcac  
acctcatcacattaa