

LOC440570-LOC440570 Gene

Reference: PTXBC126340

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	LOC440570
GENE ID	440570
ORIGIN SPECIES	Human
ACCESSIONS	BC126340
VECTOR	pDONR223
NCBI SYMBOL	LOC440570
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaatcacgaattgatgaccagttggctcttgatgtgaggaataatccagagtcagagggaactctgaagtttgcgggagcaaacggaagggtggcgttgccatgcctaaagatgggaaatggcaggtgtcacaggttcaggggaaggtcggagaccagctgaggg
ccccggagccttctggaaagagttcccatccagcccgctcgggttccgcatccgtcttattcctatgacgttgagggtgctggcgtctgggtcctttatgatgcagagggtgccccgctcaccggggcgccctccgcgtcccgctcctctggcaacctggtgcgggctccgga
tctggcgaccacgacccgctggtcacttgcctccctcgaaggcgcatctctagtcagtggtgagctgcggccgggtcgtgtaactcgtccaggactcgggactcgtggccttgggtgccctgcggagccctcgggtgtgcctgcaggtcttttttga