

ANKRD45-ankyrin repeat domain 45 Gene

Reference: PTXBC126353

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT117; ankyrin repeat domain-containing protein 45; cancer/testis antigen 117; ankyrin repeat domain 45		

For research use only.

GENE DESCRIPTION	ankyrin repeat domain 45
GENE ID	339416
ORIGIN SPECIES	Human
ACCESSIONS	BC126353
VECTOR	pDONR223
NCBI SYMBOL	ANKRD45
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtcagaaggacccacagagtcagagagtcagaattttctcacagcaagaaggagaaatgaagaagaagcccaagaaccagaggaaacagccctaagaacccctttacaacctgctctcacaggggatgtagagggtttgcagaagatattgaggatcctgagaat
cctcatcatgaacaggccatgcagcttcttagaagaagacatcgtgggagaaattgtgtatgcagcttgcagctgggcaaaagtgcagtgattagagctttggcaaaatattgtgtgaatctgaatgaaaaaccacagagggtacacactcttactgtgctgcagcctgggtcg
tttgaaactttgaaagcactggtagaactggatgtgatatagaagcttgaactccgggaagaagggtcgcagatgttgcctgctagatattctcagactgagtggtgtgaattcctggactgggcagatgcaaggctgactctgaaaaatatattgcaaaagtctctttagctgttactgaca
cagaaaagggtatcagggaactccttaagggaagacaagaatacctcagtcgcatgccgtgcataaaatgagtggttgaaacccatacagaagcttcattaatgagctttttgagcagagacaacaactggaagatattgactcctatctcacaaaaatgactacacatgtcaag
tgaaaagtgcataatctgtaacaagccatgaccagaagagaagtgcaagatgatacctccaactga