

FLJ40235-hypothetical protein FLJ40235 Gene

Reference: PTXBC126357

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C19orf75; SIGLEC23P; SIGLECP7; SIGLEC family-like protein 1; sialic acid binding Ig-like lectin, pseudogene 7; SIGLEC family like 1		

For research use only.

GENE DESCRIPTION	hypothetical protein FLJ40235
GENE ID	284369
ORIGIN SPECIES	Human
ACCESSIONS	BC126357
VECTOR	pDONR223
NCBI SYMBOL	FLJ40235
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttccactgctacagctgggtacctgctaagctgctcaactctcttgccttggagaagacactgcagtcagctgttcttccatgggattccacacctctgtgcagtggtgatgggaggagtcctccgtgggtgtggatggcatggatggcagcctcaagtgaattccaccatgc
ttggccctgggtaacagcaccatcagcctaactgaagagccagaatgggcatgagacttctctgtgagggaagaacaaaacggaacccatgctttgagcctcactgatgtcaagaagagttcttggctgccaggcctctgtgaaaggctgacccaggtgtatctatgcg
ggaattgtaattgcgctgcttctctcctcctcctctcatagtgaaacatcagaagaagcagggcgaagaagctgcagcgatcagagcaaaaagagctctaaagtcagagcaagccaagaacttgagatgtcttgaagcctgaggaaccaggaaaaccgtatgcgcacat
ttctgagagccggatattggaagcaagataaacgagccagctaa