

DEFB123-defensin, beta 123 Gene

Reference: PTXBC126365

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DEFB-23; DEFB23; ESC42-RELD; beta-defensin 123; beta defensin 23; specific to testis and epididymis; defensin beta 123		

For research use only.

GENE DESCRIPTION	defensin, beta 123
GENE ID	245936
ORIGIN SPECIES	Human
ACCESSIONS	BC126365
VECTOR	pDONR223
NCBI SYMBOL	DEFB123
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagctccttttgctgacttgactgtgctgctcttccagctgactccaggtggcaccacaaagatgctggaatctttatggcaaatgccgttacagatgctccaagaaggaaagagtctatgtttactgcataataataaaatgtgctgctgaagccaagtaccagccaaaagaaagggtggtggccattttaa