

H2BFS-H2B histone family, member S Gene

Reference: PTXBC126369

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2bfs; 2610022J01Rik; R74621; histone H2B type 1-C/E/G; H2B histone family, member S; histone 1, H2bc; histone cluster 1, H2bc		

For research use only.

GENE DESCRIPTION	H2B histone family, member S
GENE ID	54145
ORIGIN SPECIES	Human
ACCESSIONS	BC126369
VECTOR	pDONR223
NCBI SYMBOL	H2BFS
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccggagccagcgaagtccgctcccgccccaagaaggctcgaagaagccgtgactaaggcgcgagaagaaggcggcaagaagcgcaagcgagccgcaaggagagctactccgtatagctgtacaaggtgctgaagcaggtccacctgacaccggcatctccttaag
gccatgggaatcatgaactccttcgtcaacgacatcttcgaacgcatcgcaggtgaggtctccgcctgcgcattacaacaagcctcgaccatcacctccaggagatccagacggcgtgcgctgctgctgccggggagttggccaagcacgccgtgtccgagggcacaagg
ctgtcaccaagtlacaccagcgctaagtaa