

## TAL2-T-cell acute lymphocytic leukemia 2 Gene

Reference: PTXBC126373

|                |                                                                                                                                                                     |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                 | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | T-cell acute lymphocytic leukemia protein 2; T-cell acute lymphocytic leukemia 2; TAL-2; class A basic helix-loop-helix protein 19; TAL bHLH transcription factor 2 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | T-cell acute lymphocytic leukemia 2                                                               |
| <b>GENE ID</b>          | 6887                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC126373                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | TAL2                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgaccaggaagatcttcacaaataccagggagcgggtggaggcagcagaatgtcaacagcgctttgccaagctgaggaagctatccccactcaccctccagacaaaagctgagcaaaatgaaagcttcgcctggcaatgaggtatatcaactcttggtaaggtctgggggag  
caaagcctgcaacaacgggagtgctgctcaggggaacattctgggctctccctcaaggacccacctgccagcctggaggacagaactctgcttgagaactaccaggttccttcacctggtccaagccacacattccttag