

ACOT6-acyl-CoA thioesterase 6 Gene

Reference: PTXBC126378

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf42; c14_5530; acyl-CoA thioesterase 6		
For research use only.			

GENE DESCRIPTION	acyl-CoA thioesterase 6
GENE ID	641372
ORIGIN SPECIES	Human
ACCESSIONS	BC126378
VECTOR	pDONR223
NCBI SYMBOL	ACOT6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgcagcatccaaagtgaaaggtcctagattgcgtctctggattttccaaaggaggtgacctgtgtctctcaatggctctcttgaaggcagcagccactgtacttatcaatgcctgtgtgagccaacacagtagctcctctacattacaaggatatgattattcctaaactgtcgatg
atctaggaaaaagtaaaatcactaagtcaggattttctacttttatggacacttggagcaatccactggagggaacacacacacaaagtctgttccattggaaggcgcaggtgcctctctgtttattgttggcatggatgatcaaaagctggaaggtgaattctatgctcagatagcctctga
aaggctacaagctcatgggaaagaaagacccagataatctgttaccagaaactgtgacctgtattgacccacctattttctcctcttagagcttctgtgcacgctgttttgggtgaggcaataattctatggaggtgagccaaaggctcactcaaaaggcacaggtagatgcctggcagcaaa
ttcaaaactttctcataaacatctcaatggtataaaatctgtcaagcacagcaaaataaa