



Reference: PTXBC126392

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	glycogen phosphorylase, muscle form; myophosphorylase; phosphorylase, glycogen, muscle; glycogen phosphorylase, muscle associated		

For research use only.

GENE DESCRIPTION	phosphorylase, glycogen, muscle
GENE ID	5837
ORIGIN SPECIES	Human
ACCESSIONS	BC126392
VECTOR	pDONR223
NCBI SYMBOL	PYGM
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]