

C2orf76-chromosome 2 open reading frame 76 Gene

Reference: PTXBC126395

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0538 protein C2orf76; AIM29; chromosome 2 open reading frame 76		
For research use only.			

GENE DESCRIPTION	chromosome 2 open reading frame 76
GENE ID	130355
ORIGIN SPECIES	Human
ACCESSIONS	BC126395
VECTOR	pDONR223
NCBI SYMBOL	C2orf76
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctctggagaagtgacatcacagtgcctcatcgttccttgaacatcgcaatttcaaacctgtatgtatcacggagtgaatttgaccaaacgttaaaggaaattatcgtatttctaaagcaagatatccctttaaggaccaacgtgccaccaccatcagaaattataaatatgatgca
ctaaagattatcatcaagcacataaatcaagacaaatgaacttgtgtgagtttgaagatgacgaaagactcctgtgaagaagacagcactctgaaagcagctggaatcgccagtgaactgaaattgcattctctgtgaagaagattataagaactacaagctaataccattcatcc
tgggtga