



Reference: PTXBC126396

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FOXJ2 forkhead factor; forkhead box protein J2; fork head homologous X; forkhead box J2		

For research use only.

GENE DESCRIPTION	forkhead box J2
GENE ID	55810
ORIGIN SPECIES	Human
ACCESSIONS	BC126396
VECTOR	pDONR223
NCBI SYMBOL	FOXJ2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1