



Reference: PTXBC126413

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C10orf71; chromosome 10 open reading frame 71		
For research use only.			

GENE DESCRIPTION	chromosome 10 open reading frame 71
GENE ID	118461
ORIGIN SPECIES	Human
ACCESSIONS	BC126413
VECTOR	pDONR223
NCBI SYMBOL	C10orf71
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]