

## SP3-Sp3 transcription factor Gene

Reference: PTXBC126414

|                |                                                                                                                      |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                  | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                 | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | Sp3 transcription factor; GC-binding transcription factor Sp3; transcription factor Sp3; SPR2; specificity protein 3 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | Sp3 transcription factor                                                                         |
| <b>GENE ID</b>          | 6670                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC126414                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | SP3                                                                                              |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgaccgctccgaaaagcccgtgaacaagaggaaatggctgacctggacgtggatagcggcgccggcggtggcgccggcgccgacggcgagtagctgcagcagcagcaacagcagcaggaacgggtgcgggtgcggcgccgagcggcgccgagcactgacgctga  
ccgctcgtctctgctggccgtacctgcagcaagatagggccgcatcgccggcgacgagcagggaggcgccggcgccgagcggcgccggcgccgagcaggtgattggcttgcagagtaggagagcaccacacgagtgagggtttgtcagccaca  
cctacaactataaagatgaagctggtaatactagtcagatccaagtcgctactcaagtgggcagtagtcttccctcagaattggcagaatcaacaataatttccgttgaccaggatcagatcaaatggtagcgtgctcaggttcaataatcaagtataccacagatccagtc  
gcagatggtcagcaggttcaaatgggttccagggctctcagataatgggggtataaatacaagaagcagtcgaattcagatcctctgctctaatcaaaccttactgctctggaacaccttctgctaacatccagaatctcattaccacagactggtcaagtcagggttcagggtgca  
attggtggttcatcttctgctgcaaacccaagtagtgcataatgctctgctgctgccaagaataattacgtttgtaccaatcaatagtcgatctagattcttgggactctcggcgagttctcagacaatgactgcaggcattatgccgagcgacatttgataaacacaggacaagctatg  
gatagttcagacaattcagaaggactggtgagcgggtttctcctgatataatgaactaatactgatacagatttatttggccaacatcctctcattcacagttgctgttagcagtagatagtagcaggtatattacaacaaaacacaaatagctgactacatctagtgccaggttcatctcag  
atcttcagggaattatccagtcgctggttctgaaagagacacagcagacaatattcaggtttctcagcagcagccgtgtgtgacagcctcacaactcgaagagtcagcagccaacccagtcagcccaattgtgcaagggtattaccacacagacaatccatggtgtgcaagccagtggt  
caaaatatacacaacaggcttgcacaaatctcagttgcagctgaatcctggaaaccttttaattcaggcagacagcagtgacccctctggacaggttaactggcaaacgtttcaagtacaagggttcagaaactgcagaattgcaaaatcagaatactgctgcccaacaataacttgacg  
cctgttcaaacctcacttgcaggttgcggcaggtgagccttcaactcagtttagtctaagcactgtgtcagttgccaatctacaacagtttacagtgactctatagattctgctgtatagctacatccaggagagaatgctgacagtcctgcagatattagatcaaggaa  
gaagaacctgatcctgaagagtgccagctcagtggtgattctacctgaaatcaaatgacctaacacacttaagagtagcaggtgtgtagatgaagaaggggaccacaacacatcaagaaggaaagaaagacttcggagggtgactgacactgtcccaactgtaaaagaggtgtggaagaggt  
accaactctgggaaaagaagcaacacattgtcatataccaggtgtgttaagtctatgggaagacctcacatctgagagctcatctgctgtgacattctgagaaagccctttgttgaactggtgactgtgtgtaaaagatttactgaagtgatgaattacagaggcagagaagaaca  
catagcaggtgagaagaattgtttgtccagaatgttcaaacgctttatgagaagtgaccaccttgccaacatattaaacacaccagaataaaaaagggtattcactctagcagtagcaggtgctggcatctgtggaagctgcgcgagatgatacttgaattactgcaggaggacaacgcttat  
ccttgcacaaattcaacaaggctgtgttcaggataggaactgttaatacttccgccaccagcaatcaagatactcttaccacactgaaatcctttacagctgttcacagtttctggaatgagacaatggagtaa

More info: Please contact us for any question: [contact@proteogenix.fr](mailto:contact@proteogenix.fr)