



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RNF135-ring finger protein 135 Gene

Reference: PTXBC126420

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

E3 ubiquitin-protein ligase RNF135; L13; MMFD; REUL; Riplet; RIG-I E3 ubiquitin ligase; ring finger protein 135

For research use only.

GENE DESCRIPTION	ring finger protein 135
GENE ID	84282
ORIGIN SPECIES	Human
ACCESSIONS	BC126420
VECTOR	pDONR223
NCBI SYMBOL	RNF135
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr