



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

KRTAP4-7-keratin associated protein 4-7 Gene

Reference: PTXBC126421

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP4.7; KRTAP4.7; Ultrahigh sulfur keratin-associated protein 4.7; keratin associated protein 4-7		

For research use only.

GENE DESCRIPTION	keratin associated protein 4-7
GENE ID	85287
ORIGIN SPECIES	Human
ACCESSIONS	BC126421
VECTOR	pDONR223
NCBI SYMBOL	KRTAP4-7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]