

BEX1-brain expressed, X-linked 1 Gene

Reference: PTXBC126427

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein BEX1; BEX2; HBEX2; HGR74-h; brain-expressed X-linked protein 1; ovarian granulosa cell 13.0 kDa protein hGR74; brain expressed X-linked 1		

For research use only.

GENE DESCRIPTION	brain expressed, X-linked 1
GENE ID	55859
ORIGIN SPECIES	Human
ACCESSIONS	BC126427
VECTOR	pDONR223
NCBI SYMBOL	BEX1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagtccaaagagaacgagcagtaaacagtcagcatgaaaatgccaaagaaaatgaagaaaggagcaagttgctaataaaggagcccttggccctcccttggatgctggtgaatactgtgtgcctagaggaaatcgtaggcgggtccgcgttaggcagcccatcctg
cagtatagatgggatgatgatgcataggcttgagaaccacaggcaaggatgagagaagagaatggaaggattggggaggaggtgagacagctgatgaaaagctgagggaagcagttgagtcagtcctcgggcagtcagcactgacccccctcaccatgacatcatgat
gagttttgccttatgccctga