

GNRH1-gonadotropin-releasing hormone 1 (luteinizing-releasing hormone) Gene

Reference: PTXBC126437

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GNRH; GRH; LHRH; LNRH; progonadoliberin-1; GnRH-associated peptide 1; gonadotropin-releasing hormone 1 (luteinizing-releasing hormone); leuteinizing-releasing hormone; luliberin I; prolactin release-inhibiting factor; gonadotropin releasing hormone 1		

For research use only.

GENE DESCRIPTION	gonadotropin-releasing hormone 1 (luteinizing-releasing hormone)
GENE ID	2796
ORIGIN SPECIES	Human
ACCESSIONS	BC126437
VECTOR	pDONR223
NCBI SYMBOL	GNRH1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagccaattcaaaaactcctagctggccttattctactgacttggtgcgtggaagctgctccagccagcactggctctatggactgcgccttggaagaaagagagatgccgaaatttgattgattcttccaagagatagcaagagggttggtcaactggcagaacccaacgcttc
gaatgcaccacgcaccagccagcttctccctccgagacctgaaaggagctctggaagtctgattgaagaggaaactgggcagaagaagatttaa