

UTS2-urotensin 2 Gene

Reference: PTXBC126443

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRO1068; U-II; UCN2; UII; urotensin-2; prepro U-II; urotensin II; urotensin 2		

For research use only.

GENE DESCRIPTION	urotensin 2
GENE ID	10911
ORIGIN SPECIES	Human
ACCESSIONS	BC126443
VECTOR	pDONR223
NCBI SYMBOL	UTS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatttattgctcacctcatgaagacgcgccttaactccggaggagctagaaagagcttccctctacagatactgccagagatgctgggtgcagaaagagggatatttcaggaaagcagactcaagtaccaacattttaaccaagaggaaattgagaaagttcaggatttctctg
gacaagatcctaacaattttactgagtcattttggccagaatctggaaaccatacaagaacgtgagactcctgattgcttctggaaatactgtgtctga