



Reference: PTXBC126445

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC065; coiled-coil domain-containing protein 113; coiled-coil domain containing 113		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 113
GENE ID	29070
ORIGIN SPECIES	Human
ACCESSIONS	BC126445
VECTOR	pDONR223
NCBI SYMBOL	CCDC113
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]