



Reference: PTXBC126451

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PP1G; PPP1R3; protein phosphatase 1 regulatory subunit 3A; RG1; glycogen and sarcoplasmic reticulum binding subunit, skeletal muscle; glycogen-associated regulatory subunit of protein phosphatase-1; protein phosphatase 1 glycogen- associated regulatory su		

For research use only.

GENE DESCRIPTION	protein phosphatase 1, regulatory (inhibitor) subunit 3A
GENE ID	5506
ORIGIN SPECIES	Human
ACCESSIONS	BC126451
VECTOR	pDONR223
NCBI SYMBOL	PPP1R3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]