

PRR3-proline rich 3 Gene

Reference: PTXBC126457

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAT56; proline-rich protein 3; MHC class I region proline-rich protein CAT56; proline-rich polypeptide 3; proline rich 3		

For research use only.

GENE DESCRIPTION	proline rich 3
GENE ID	80742
ORIGIN SPECIES	Human
ACCESSIONS	BC126457
VECTOR	pDONR223
NCBI SYMBOL	PRR3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgaacgaagaagcagaatcatcaccagccaccgacacagcagcagcccccgtgccgagcgggaagagactggagatgaggaggatgggagtcctatcgaccaccagcctctgggcctcccccatggccaatggaaaacctggcgaccctaagtcagctcttc
acagaagtctccaggatcaaggaggaccactgattccaccactgctgagtcctccacctctctctgggtagaggcccaattcggagagggcttggccccagggtctagcccatatggtctgtgtgtgggagtcgaatgcagaacctctttccggggccaggccatgggggtccca
ccaggggaagctttcacaaagaacagagaacctcgaaggctcaaaagctggctcttcaagaatacctgccgcccaaggatgaccccaagggttatggaagacaaatccgaccgcctgtctgccgacatttgcacaaaaggggccactgtcgatagaggacctctgtgccttcta
ccaccagcgtcaatggacctctctgtga