

SFRS12IP1-SFRS12-interacting protein 1 Gene

Reference: PTXBC127085

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein SFRS12IP1; SFRS12IP1; P18SRP; protein SREK1IP1; SFRS12-interacting protein 1; p18 splicing regulatory protein; splicing regulatory protein of 18 kDa; SREK1 interacting protein 1		

For research use only.

GENE DESCRIPTION	SFRS12-interacting protein 1
GENE ID	285672
ORIGIN SPECIES	Human
ACCESSIONS	BC127085
VECTOR	pDONR223
NCBI SYMBOL	SFRS12IP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagtcgccaggttgcaacaaggacagtgtcagagcagcgtgtgtaaaaaatgtggtacccgtgtcacctgacttttgatgccgcaattttccgagtagaccctaaaaggacatagttttggtgtcagcagtagcaagtagtgagatagcgtgaagagaatgaagaactgaataaa
ttgcaggcattacaggaaaaagaataaatgaagaggaagaaagaaaaagcaaaagaaaaaatcaattgaaaaaaaaggaaaagggtcttctcatccagttccactgaagaggacactcaaaacaaaagaaacaaaataatcagaagaaagaaaaagaa
aaaaagagtaaatcaaaaaagggaacatcacaaaaaggaaaaaagaagagaaaaaggaaaaagcatttctacacctaattgttctgaattctcagaagtaa