



Reference: PTXBC127102

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSPC082; SIP-1; SIP1; SMADIP1; ZFHX1B; zinc finger E-box-binding homeobox 2; Smad-interacting protein 1; zinc finger homeobox 1b; zinc finger E-box binding homeobox 2		

For research use only.

GENE DESCRIPTION	zinc finger E-box binding homeobox 2
GENE ID	9839
ORIGIN SPECIES	Human
ACCESSIONS	BC127102
VECTOR	pDONR223
NCBI SYMBOL	ZEB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr