

RDM1-RAD52 motif 1 Gene

Reference: PTXBC127190

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAD52B; RAD52 motif-containing protein 1; RAD52 homolog B; RAD52 motif 1; RAD52 motif containing 1		

For research use only.

GENE DESCRIPTION	RAD52 motif 1
GENE ID	201299
ORIGIN SPECIES	Human
ACCESSIONS	BC127190
VECTOR	pDONR223
NCBI SYMBOL	RDM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggagttggtaccttttgcgggtccatcgagagtgacaaaaccttgcctgtgtggagctgagctccggaccacggccgaggcttgcacattctctgttcacagcattttcagtttggcctctgtattcagtcgggtcttcccaaatgctgcagtgcccatcctgtttctatgc
cgtcattaagtgtttatctgcaagggtgcccacagagcccaaaagcatgcgaccggaaagcagctttttcagaaatctccagtcagggtctgttggcaccagacataaggcagttcaacatcaagccctgacagttccaaatgccaagaactggcgaattactactttgtttca
atgggtgttccaaaaggatcatcaagcttcaggagcttctgaccttgaagaaagggaatgaagatagcatggtgccacttccgaagcaaacctgaagttctctgtgtctttagaagtggtgttgcctcctgtgattgcaggagtcctggcattggcttggaggagcctatggataa
ggtggaggagagaccattatcattccttgaaggaagacagcccaagcttgcattcagaaggcttgcagatgcattccagaaactgttgattgtttctagaaagtggtgaaatagctgtggagtacagaccagtgagacatcgtaggtgtcagatgcgaagaagaactacac
ggtttaattcaagtccttgcctcctggaagcagtaggccaagagggaagggtatctctcggattcagcttgaggagggaagagtcaggctgccagaacttgactag