

HIST1H3C-histone cluster 1, H3c Gene

Reference: PTXBC127610

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H3.1; H3/c; H3FC; histone H3.1; H3 histone family, member C; histone 1, H3c; histone H3/c; histone cluster 1, H3c; histone cluster 1 H3 family member c		

For research use only.

GENE DESCRIPTION	histone cluster 1, H3c
GENE ID	8352
ORIGIN SPECIES	Human
ACCESSIONS	BC127610
VECTOR	pDONR223
NCBI SYMBOL	HIST1H3C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcgtacgaagcaaacagctcgaagtctaccggcgcaagctcgcgaagcagcttgcactaaagcagccgtaagagcgtccggccaccgggtggcgtgaagaaacctatcgctaccgcccgggcaccgtggccttgcgcgaaatccgtcgtaccagaagtcacc
gagctgctgatccggaagctgccgtccagcgtggtgcgagaaatcgccaggacttcaaaacgacctgcgtttccagagctctgcggtgatggcgtgcaggagcttgtgaggcctacctgggtggactcttcgaagacaccaatctgtgcgtattcagcgtaaacgcgtcacc
atcatgcccaagatattccagctggcagctcgcacatccgtggggaaagggcataa