



Reference: PTXBC127622

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KMLCK; MLCK; MLCK2; skMLCK; myosin light chain kinase 2, skeletal/cardiac muscle; myosin light chain kinase 2, skeletal muscle; myosin light chain kinase 2		

For research use only.

GENE DESCRIPTION	myosin light chain kinase 2
GENE ID	85366
ORIGIN SPECIES	Human
ACCESSIONS	BC127622
VECTOR	pDONR223
NCBI SYMBOL	MYLK2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr