

SPATA1-spermatogenesis associated 1 Gene

Reference: PTXBC127634

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	4921536I21Rik; AV266943; SP-2; spermatogenesis-associated protein 1; spermatogenesis associated 1		

For research use only.

GENE DESCRIPTION	spermatogenesis associated 1
GENE ID	64173
ORIGIN SPECIES	Human
ACCESSIONS	BC127634
VECTOR	pDONR223
NCBI SYMBOL	SPATA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagaagataaaaaagaacaaattccaattgataattatccattcagacattagtgatgtcactcaatccaagtcgaccttcctcatcagagttggtggaactcatgtttttatgtccctgaaggatcatggaactataagctaaataccattcaacagaagtgttaacaaattcat
ttcagctggatttctaagagtatctcctcaacttctacgagccctgagggagcgtcttggtgagttcctgggtgaagatgctattgcagaaaaattttattctgaaatgcattggaataatttagctgtggtggaaggaagcaagaatcagaactgaaactcaaatcattgtctctccat
gctcttcaaccagaattatattgcttctgtaatggacatttaggaaatgtttatcaccatcaacagttatttttagatgagcggcagactaataatggtttaatgagcctgatggaacaatccacagaccaattagtgtaactttgtcaaggaggaaactggaagatcccagttgttagaaa
acactttgaaagagcttctaacaagaatcaggagaagcgtgggggaaagccactgcagaaaaagccaaattgcaaaaaatcaaatggaattctgagttgccagatcattggaagattcaataatgattgctttggcactaaaaaagtgtctttgggaaatgaagatgatacagc
tatcagtagaagacaggacaatcagacagctgaaaaagagtcacacccctaccagatcaccttcacttctgtcaacctgttctttcaggaataa