

LOC285141-similar to CG14853-PB Gene

Reference: PTXBC127676

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to CG14853-PB
GENE ID	285141
ORIGIN SPECIES	Human
ACCESSIONS	BC127676
VECTOR	pDONR223
NCBI SYMBOL	LOC285141
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagactgtaaatgaaccagaacaggtgaagtgaagcagaagatgcagtcattgtaaagcaggagaaaaataatgaattgccttcaggatattgatgataaattgcagaatcagcagaggatgatggatgaagatgataccaatgatgaagatgatgaagatagtaaccctaaaaa
gaatactcaggccccactagagttaatggcagaattctgagagcagaatggcccgagagatcagctggcaaaaaattatgcagatgactctaatctatgaaccagaaaatcctgaggccaaggagttttcacacttattgaagaatgttctgatggagaaactcagaatcatgag
caagacggtgaaaacagtgatgaagacagcagtggtgagagtaaggagaaagcagtgaggagctgagtgacgagagctctgacgaagggtgaagatggatcatga