

MGC3196-hypothetical protein MGC3196 Gene

Reference: PTXBC127679

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 223		
For research use only.			

GENE DESCRIPTION	hypothetical protein MGC3196
GENE ID	79064
ORIGIN SPECIES	Human
ACCESSIONS	BC127679
VECTOR	pDONR223
NCBI SYMBOL	MGC3196
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgcccttgaggcgatggccacggggctgtagccgtgctgcccctgctcactgcccgcctgcaaggcagcagctgcaacgggatgtgtgctctttgagcatgatggggccgctttcaccatcctcgggctgttctgcgcgggaggggcgtctctggc
ttccatggctgtgagccgtgtccggccccgggtccgggtgcagcctctgtagcgagggtcccaaatcgtggcccttgacgtgcgtccgctctggcgctacggtctggccgtggctgcggcgccatcgagccctcgactcgggtgctgtctctctctccggtctgt
gcgctcagtggtgcttcgagctggaggcagcaggtgacctcaccactcatgcccccttggcttggggcccaattcacagttccttgaagcaggtatcttgcagtgcccaccgggggtgaagtcctgcatgctaccttgaaagtcaaggccgacgcttctattcctcttgacaaa
actggacacttccctaacacaaaactctttgacaatactgtgggtgcctaccggagcttgtga