

DAPL1-death associated protein-like 1 Gene

Reference: PTXBC127682

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	death-associated protein-like 1; early epithelial differentiation-associated protein; death associated protein like 1		

For research use only.

GENE DESCRIPTION	death associated protein-like 1
GENE ID	92196
ORIGIN SPECIES	Human
ACCESSIONS	BC127682
VECTOR	pDONR223
NCBI SYMBOL	DAPL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaatgaagtgcagacctgctctccctcggaagggggacatcctcgcagtaaaagctggaggaatgagaatttccaaaaacaagaattggcacctggaaagacataccaaaaacaggattcgagaaaaaagtgccattgcaaatgttgccaaatacagacacc
ggatgccctgaatgacgcactggagaagctcaactataaattccagcaacagtgacatggcacatcaaaaacccacacctgctctggaaaagggtgtccactgaaaaggatctacattattcagcagcctcgaaaatgttaa