

SHISA3-shisa homolog 3 (Xenopus laevis) Gene

Reference: PTXBC127690

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	hShisa3; protein shisa-3 homolog; shisa homolog 3; shisa family member 3		
For research use only.			

GENE DESCRIPTION	shisa homolog 3 (Xenopus laevis)
GENE ID	152573
ORIGIN SPECIES	Human
ACCESSIONS	BC127690
VECTOR	pDONR223
NCBI SYMBOL	SHISA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggcactgctggcgtttgcttctccttgctgtctgctggggccggggcgccagcagtcggagagtagtgcacggctgggtggagctgcagggaactaccagagggttcagtgccagaggacttcacagctggagctaccatctgctggctcct
gcgcgtccgctactgttgcggccggcgagcggaggtggagcaggcggtgacacagaccgcccgaactggagcaccagggcatcactgcgcagcctgtctacgtccctttctcgtcgcgtccatcttcattgcgttcacatcctgggctctgtatggctatttattgtt
gcacctgttggagaccgaaggagcctcgcagcagccaatccgcttctcactccgcagctatcagacagagacctgcccagatcctgacctccacagccccagggcaccctccggcagtcacagccagagctccagctccacagggcctccatccgaggttctccttt
gccagggtgagccggctgctgtggtgcctcaccgccccgcatacaccagcactcaatccactggctcagccatctggttctctgtgtacccagctatttgcgttaccctccagcaggagccccactgctgggaagagctgtccagacttcagttcagttga