

TMEM220-transmembrane protein 220 Gene

Reference: PTXBC127705

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 220		
For research use only.			

GENE DESCRIPTION	transmembrane protein 220
GENE ID	388335
ORIGIN SPECIES	Human
ACCESSIONS	BC127705
VECTOR	pDONR223
NCBI SYMBOL	TMEM220
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccagcgctgtggcgggctgcaacggactcggcgcccttctcgcgctggcgcccttggtgcaggtaaatgacccagatgcagaggtgtgggtggtgtgtacacaatccctgcagtactgacctgcttggacttaacctgaagtcacaggtaatgttattggaaaag
tatctctgcaatacacatactctttgtacgggtgtggcgtgtggcgtgctcctacctctgcatcgtacacaacagaaacatcttacatgaggaagaaggcaggagcgtctctggtgattattacagcatggattatcctgtgccacagttcctcaagaatccagttgtgggaagaatc
aattggctattgccattgtaatacacatttcccaattatctcatgggtctacatatataaacaaggaaatgcggtcctcttgccaactcactgcaagacagtaatttaa