

BCL2L15-BCL2-like 15 Gene

Reference: PTXBC127719

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf178; bcl-2-like protein 15; bcl2-L-15; pro-apoptotic Bcl-2 protein; BCL2 like 15		
For research use only.			

GENE DESCRIPTION	BCL2-like 15
GENE ID	440603
ORIGIN SPECIES	Human
ACCESSIONS	BC127719
VECTOR	pDONR223
NCBI SYMBOL	BCL2L15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagagctccaaacttttgaggacaacgaatgcattgtgaacctctactcaggactcttgagcccaacattgcaggtgccagccgggaacatgctgttagatgaagtagattcaggagagcctgttcttttgatggctatcattgctgctgccttcggatgttggtgacc
agttcaacggagaattggaagcttctgccccaaacgtcattgcagaaccattaaaggacagacagagctatactccagaacactgtggaatctctcagaagacctggtgtgctcaggattccagcttagcttatgagagagctttctggcagtgctcagtgaaactcttgatcatgct
cacattgctcctgaagtagtgggacaggtggctatccccatgacgggtatgatcaatgggaaccaagccatccgggagttcatcagggccaggagggtgggaaatctggagagctga