

**C1orf130-chromosome 1 open reading frame 130 Gene**

Reference: PTXBC127785

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf130; MP11; noncompact myelin-associated protein; myelin protein of 11 kDa; non-compact myelin associated protein		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	chromosome 1 open reading frame 130
<b>GENE ID</b>	400746
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC127785
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	C1orf130
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

**ORF (Only for reference, product provided without STOP codon )**

atgaccacagccaccctctggggataccacctctctcactgaacatgaccaccaggggagaagactctctgataagattctggagccattgttctgcccgttgggtggtgcatcatcatctcaccgtggttctgctcctgctgaagatgtacaacagaaaaatgaggacgagcc  
gggaactagagcccaagggcccaagccaaccgccctctgcccgtggcccaaacagcaacggcagccaacccagcaactgtgacctcagctctgtgacgtccaggtggagacgcatga