

TMEM72-transmembrane protein 72 Gene

Reference: PTXBC127790

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf127; KSP37; transmembrane protein 72; kidney-specific secretory protein of 37 kDa		

For research use only.

GENE DESCRIPTION	transmembrane protein 72
GENE ID	643236
ORIGIN SPECIES	Human
ACCESSIONS	BC127790
VECTOR	pDONR223
NCBI SYMBOL	TMEM72
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagctccaggtgtctgactgggctggaatacacctgcggctcctgggcatcaccactgctgcagtggtgacgcgtgggcaactgagacctctccaggccagttcaaaagcctggctttctatctgctgtttacaggagccgctgtctccatatgtgaaggggctactttgtgg
ctcagctgctggccatctgtctccaggtgtcaaccagggtccctggcagacagagtaaggagaaagccactggctgggctgcttccagaagttcctggcctaccctgctgctgctgggctgttctctccaccggctcgtgtgtgacgtgacctccaggctccatgctcatc
accggcctggcctacttctctgagcaagcgggaagaagagaaagctgccccgaggtgctggcctccccagagcagtagacagacccctctagcagcgtgtgagcaccaccggctctggggacacagagcaaacctacactttccatggggccctcaaggagggccagctc
ccttttaccacatgaagagtagctgaaggggactaagaagccagtgccctccagcccccaacacctgatggagctgagcctggagccagccgactccctggccaagaagaagcaggtgcactttgaagacaacttggctccgcatagtcctccctcctcgtgaaggtctgga
tggggacagtgagccagaggagaccacctctgacacgacacccatcctccctcccccagggcccccactctctctgctcatctctacagccaccgctgtctga