

KIAA0574-KIAA0574 protein Gene

Reference: PTXBC127817

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TMEM228; protein FAM189A1; transmembrane protein 228; family with sequence similarity 189 member A1		

For research use only.

GENE DESCRIPTION	KIAA0574 protein
GENE ID	23359
ORIGIN SPECIES	Human
ACCESSIONS	BC127817
VECTOR	pDONR223
NCBI SYMBOL	KIAA0574
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtgtatgcagatgtctctccgatgtctcgcagatgtctctccgcagaggtcacatctgccaacccacatgctgactctcagggcacgttctccacagacccctggatttcagagttcatccccactcccgctccaccctattacccccagagtacacctgcacac
cggtccactgagggccagaggggctccacctggacttgcctcgtctccattcggcactttgtatgacgtggccatcaatagccccggcctgctctatctctgagctccccctccgtacgaggggtggtggccagccccctccagccaggttacaagtataggtcagcaggtggc
cgagtcagctccgggacccaacaccagtgctggttcagcactccagttaccagctgacagcacaagcctctgtgtccgagggcactgctacgaggttccagccatccccgatggcctgtggcgccccagcaccctccgaactgcccctccccggccgctgtctc
cagaggatcctggcctgacaggtgcagccaggtccgggcatgtgtcccgctccaccagcgacccacactgtgacacatagcatggcagcgcatgcctcttcacagggcgtgtgttagccaggacctggaagcgggactgtctgaggtgtgcctggagcgttccat
gtctcgtctgccacggctgctgctgggccagcttcaccagcggggacccagatactggaactgaccaacggccaacaccagagccctccagcgacctccaagagcggccacgctctttag