

C17orf89-chromosome 17 open reading frame 89 Gene

Reference: PTXBC127837

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C17orf89; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 8; NADH:ubiquinone oxidoreductase complex assembly factor 8		

For research use only.

GENE DESCRIPTION	chromosome 17 open reading frame 89
GENE ID	284184
ORIGIN SPECIES	Human
ACCESSIONS	BC127837
VECTOR	pDONR223
NCBI SYMBOL	C17orf89
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcggctaacggagcgggtgtgtggcgccgtgcgaagccgctccgccttccccgagcggctggccgctgcggggccgagggccgctacggcaggtgcgtgcaggcctccacggccccggcgccgctgagtaaggactctgcgcggggagttcgagccct
gcggagctgcttcgccgctgcggccaagaagacgctggaggaggctgtag