

C22orf41-chromosome 22 open reading frame 41 Gene

Reference: PTXBC127859

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C22orf41; THEG2; synaptonemal complex central element protein 3; testis highly expressed gene 2 protein; testis highly expressed protein 2		

For research use only.

GENE DESCRIPTION	chromosome 22 open reading frame 41
GENE ID	644186
ORIGIN SPECIES	Human
ACCESSIONS	BC127859
VECTOR	pDONR223
NCBI SYMBOL	C22orf41
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgatgctgacccctgaggaaagaactatgacaacatgctgaaaatgctgacagatctgaataaggacttgaaaagctattagaagagatggagaaatctcagtcaggcgacctggatggcctatgacatggtggtgatgcgaccaaaccctacgtggccgagtcctatgcgtc
ggctggaggatgccttcgtaactgcaaggaggagatggagaagaactggcaagagctgctgcatgagaccaagcaaaaggctgtag