

CLRN2-clarin 2 Gene

Reference: PTXBC127863

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	clarin-2; clarin 2		
For research use only.			

GENE DESCRIPTION	clarin 2
GENE ID	645104
ORIGIN SPECIES	Human
ACCESSIONS	BC127863
VECTOR	pDONR223
NCBI SYMBOL	CLRN2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcctggatgttcaaaaagcgtggtatggcgtggcgttttactcagcttctctcttctcctgatcatcgttgcctggtagtgtcccaactggctgagtggaataatcctttgcagactggagtggatctgtcaacgccacagacagagagctggtaagttcattggggacatttac
tacgggctctccgaggtgttaagggtcggcagtggtggcgttggggccgccaatccaattcacgattctccacacctggtgaaggagctcaacgcagccttcattgtgatgtctgtctctctctctggccttggccctggtctgtcagcatggccttggcattcttaacatgat
ccagggtccgtaccggcagtcagcgttctggggcattctgctatggaatgtctggtcagcggcgtctgtgctgttagccatcggcagcttctgtgctgctggtgaaatttcacgacctgacggaacgaatcgccaacttcaggagagctcttcagtttgggtgggaagaaca
gtatgaagagtcgttttgatctgctggtggcagcgttcggccatgctgcaaaactgtgctgtggtggcagtcagtaaatccctccctgagattaagacaaaatgaagagggccacggtcacagctgaggatattctgtattaa