

WDR38-WD repeat domain 38 Gene

Reference: PTXBC127949

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WD repeat-containing protein 38; WD repeat domain 38		
For research use only.			

GENE DESCRIPTION	WD repeat domain 38
GENE ID	401551
ORIGIN SPECIES	Human
ACCESSIONS	BC127949
VECTOR	pDONR223
NCBI SYMBOL	WDR38
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacagcggggtccggccacgtggccgtgcccagagtgaaattcttccggccagcggcgagggtcaactcttctgcttctccctgatggccagatgctgctacaggctcagaggatggctgctgtatggctgggagaccggagtgggcagctgctgtggagctgg
gtggccacacagggcccccgtgaagtctgcccgttctccctgatggccacctcttccgagcgccctctgtgactgcaactgctgctgggatgtggcgaagagcgaagtgtctgctgggtcctgaagggtcaccacggagtggtgagacgggtcagctcagccctgactcagacagc
tgccatcagggtgctgggacaaagcgggtgatgctctgggatgtgcagctccgcccagatgtgcgcctcttagttgggcaccgtgactccatccagagcagcgaactctcaccacgggtgaactgctggccaccggctcctgggactccaccgtacgcatctgggacctgaggatgtg
acccacagctctccaccaggcgtagaggggacacagtggaacatcagctgctgtgctattcagcatccggcctctggcctcggctcctgggacagacacatcctggaagcccaaccacagcagctgcttatcgaactgaaggccatgtcacctgggtgaagagcat
agccttctctccgacgagctgtggctggccagcggcgtattccgcagtgctcaagtctgggactgcaacacaggaagtgccttgagacctgaaggagctcctggtgtggccacacctgtgcttaccaccagatgggaaatcttagtctgtgagctgcccagactaga
cgtcaaatatcccgcagctccaatcaccaggacacctcaaacctaa