

H2AFB1-H2A histone family, member B1 Gene

Reference: PTXBC128034

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2A.Bbd; histone H2A-Bbd type 1; H2A Barr body-deficient; H2A histone family member B1		

For research use only.

GENE DESCRIPTION	H2A histone family, member B1
GENE ID	474382
ORIGIN SPECIES	Human
ACCESSIONS	BC128034
VECTOR	pDONR223
NCBI SYMBOL	H2AFB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgaggaggagagacgccgagggtcctccgtgctggcggccggggcgacgtctctgcaccgtccgagcggagcttctgtttcagtgagccaggtggagcgcagtcacggaggccactacgtcagcgcctgagtcgcacggcgccggtctacctcgtcggt
tattgagtacgtgacggccaaggtcccgagctggcgggcaacgaggccagaacagcggagagcgggaacatcactccctgctgctggacatggtgttcacaacgacaggctactgagcacccctttcaacagaccacctctctaaagtggccctggcgaggactag