

PLXNB1-plexin B1 Gene

Reference: PTXBC128094

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PLEXIN-B1; PLXN5; SEP; plexin 5; semaphorin receptor SEP; plexin B1		

For research use only.

GENE DESCRIPTION	plexin B1
GENE ID	5364
ORIGIN SPECIES	Human
ACCESSIONS	BC128094
VECTOR	pDONR223
NCBI SYMBOL	PLXNB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccattccggcacaagcctgggagtggtgtctccgtggagggggagaaacctggacctgcaatgtccaaggaggagtggtggctatgataggggatggccctgtgtgtgaagacgtgacgcggcaccacctgtactgcgagccccctggagcagcccctgccacggcac
catgccctccgagaggcacctgactcttgcctgagttcacggtcagtgaggcaggtccctggccgggctgggcattga