

FLJ40434-hypothetical FLJ40434 Gene

Reference: PTXBC128116

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical FLJ40434
GENE ID	163742
ORIGIN SPECIES	Human
ACCESSIONS	BC128116
VECTOR	pDONR223
NCBI SYMBOL	FLJ40434
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatgaagttcgcttcgagcgcacggtagaggacctctacaggtatgccgtgcccaagcccagagccagtgccacaaaggccaagcagctggggcatcaccttcgtggcgcttctgcgctgtcgtgtccatcccgccgactccgtgggtgtctgtgtga