



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

GPR103-G protein-coupled receptor 103 Gene

Reference: PTXBC128133

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GPR103; AQ27; SP9155; pyroglutamylated RFamide peptide receptor; G-protein coupled receptor 103; QRFP receptor; orexigenic neuropeptide QRFP receptor; peptide P518 receptor		

For research use only.

GENE DESCRIPTION	G protein-coupled receptor 103
GENE ID	84109
ORIGIN SPECIES	Human
ACCESSIONS	BC128133
VECTOR	pDONR223
NCBI SYMBOL	GPR103
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr