



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

LUZP4-leucine zipper protein 4 Gene

Reference: PTXBC128134

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

CT-28; CT-8; CT28; HOM-TES-85; leucine zipper protein 4; HOM-TES-85 tumor antigen; cancer/testis antigen 28; tumor antigen HOM-TES-85

For research use only.

GENE DESCRIPTION	leucine zipper protein 4
GENE ID	51213
ORIGIN SPECIES	Human
ACCESSIONS	BC128134
VECTOR	pDONR223
NCBI SYMBOL	LUZP4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcttcctgtttcggagactaacgcttcttctgaaaagtgcgcgccaatacatccagctcggaaaaaggttaactcttagatattgtctctagacgattataactataaagagttagaagggacaaatgctgaagaagaaaaaataaaagacagaaccatagtaaaaagggaatcgcttcaa
gacgacaatcaacagatcagatcgcctatcgccatcgaggagcgtactcaagatcaggagaagactctgaggagaagaatctagataaaaacatcccaaaaactcttcattcgaatgctggacaacaccccttaaatggcgccgctcttaattgacagggaagtgctgacaaattaga
ggcccaacgacagagaactcaacggcctctagagggggaacagcagatcagagagaagaatctcgacaatcagaagaatccggcgccaacacagaagaataatctattctgctgagctccgcaaacactcttagagatctctcttcagtcagacagatctcaaggcgactaaa
gagacatcatcccaataatgagagatctcatggccaatacaagatctcatgttcaactctgagagatctcatgtccactcagagagatctctgctcactcagagagatctctgtgctactcagagagatctcatgtcactcaagaagatctctgtagccaggagatctctgacactca
gagtgatctcatgacctcagagagatctcatcgaactcagaaagatctcatagccactcagagagatctcatgcccactcagagagatctctgtgcccactcagagagatctcataatcagtcaggggagatctcatggccaactcagaagaacatcagaga
gatactcaacagctaaaaatacaataactataa