

SPINK5L3-serine PI Kazal type 5-like 3 Gene

Reference: PTXBC128164

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPINK5L3; HBVDNAPT1; HESPINTOR; LiESP6; serine protease inhibitor Kazal-type 13; Hepatitis B Virus (HBV) DNA polymerase transactivated protein 1; hepatitis B virus DNA polymerase transactivated human serine protease inhibitor; hepatitis B virus DNA polym		

For research use only.

GENE DESCRIPTION	serine PI Kazal type 5-like 3
GENE ID	153218
ORIGIN SPECIES	Human
ACCESSIONS	BC128164
VECTOR	pDONR223
NCBI SYMBOL	SPINK5L3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgctttccccacaagattatatttcttggtatgctctactttgacacatgtggctttctcaggaaatttcaataaacgtgacttcactagtggtcctaagccccgatgtaaaatgtatatccactggaccctgattacaatgcagactgccccaatgtgacagcacctgtttgtgcctcaaa
tggccacactttccagaatgagtggtttctttgttgaacagagggaatttcattatctataaaatttgaaaaatggaatgtgattaa