

IRGM-immunity-related GTPase family, M Gene

Reference: PTXBC128168

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Irgm; Ifggd3; Ifi1; Iipg3; Iipg3; LRG-47; immunity-related GTPase family M protein 1; LPS-stimulated RAW 264.7 macrophage protein 47; OTTMUSPWKG00059374; OTTMUSWSBG00059273; immunity related GTPase m1; immunity-related GTPase family, M; interferon-gamma-i		

For research use only.

GENE DESCRIPTION	immunity-related GTPase family, M
GENE ID	345611
ORIGIN SPECIES	Human
ACCESSIONS	BC128168
VECTOR	pDONR223
NCBI SYMBOL	IRGM
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatgttgagaaagcctcagcagatgggaactgccagagtgatcttaacatcaaggagactctgaagatagctccaggacaccagttaacatcactatggcaggggactctgcaatgggatgtccacctcatcagtgcccttgaacacacaggacatgaggtaaggcctcac
ctctactgagctgtgaaaagctaccacaagatgtgctctatttcttcccacttttcaaatgtgggtgtgggacctgctggcacagggtctgccaccacaacctggagaactactctgatggaatgcagttcaaccggtatgactcatcatggttgcacatgcacaatcagcatgaa
tcatgtgatgcttgcacaacccgtgaggacatgggaaagaagttctacattgtctggaccaagctagacatggacctcagcacaggtgccctcccagaagtgcagctactgcagatcagagaaaatgtcctgaaaatcccagaaggagcgggtatgtaatactaa