

C21orf125-chromosome 21 open reading frame 125 Gene

Reference: PTXBC128187

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C21orf125; NCRNA00319; PRED49; long intergenic non-protein coding RNA 319		
For research use only.			

GENE DESCRIPTION	chromosome 21 open reading frame 125
GENE ID	284836
ORIGIN SPECIES	Human
ACCESSIONS	BC128187
VECTOR	pDONR223
NCBI SYMBOL	C21orf125
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgttgcctcgggacacattaggaatgaaccaggcctcagtcctggcctcctcctcaaggctgggtgctcagacgtccccctcagcccttctctagaccttcggctgaacaaacgagaacatctcacagccctgcaacagccaggagctgctctggagctccacgccggct
ggagagagaggagcaccctggcagcccgccagcaccaggctctcacagtggagccacgtctcttggcctctatttgaagctgttttccacagcgcttggacaagaatgatgtgagtgagccggcgtgtcaatagggccagcacaggactggcctcaggccacattg
gggaggggtggcgaccaggccctggctgcagccacgtgggcaaaagcatggccaggcaggaggtgggcctgtgagcttggctacagcgaccgtgggcgtgggccccgggctcctcagccctga