

TMEM190-transmembrane protein 190 Gene

Reference: PTXBC128189

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MDAC1; transmembrane protein 190		
For research use only.			

GENE DESCRIPTION	transmembrane protein 190
GENE ID	147744
ORIGIN SPECIES	Human
ACCESSIONS	BC128189
VECTOR	pDONR223
NCBI SYMBOL	TMEM190
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgggtgtggatccagcgtggcctgtctctgtgtcagggctcggcagacggaatggaatccaggatcttctacccatggagctgtgagggtgacatatgggacggggagagctgtggggccaggcggccatcgaagcccaacctctgcctgcgtctcgggtg
ctgtaccgcaatgggtgtctaccaccagcgtccagacgaaacgtcggagggaagcacatgtggcgctgtgttgacgtgcagggcctctctctctgagctgcagcatctgttctgtgtggccaagcggcgacgtgtcatatgccgggttctctggcggtccgt
gtgacatgtccaaagtcgtctgtctcctcaagcaccgagggaccaagaagacgccgtccacgggcagcgtgccagtcgccctgtccaaagatgccaggatgtgaggaggaccaggggggaaggacggaggaggtgaggagacagagggcgagggaaggagg
attag